4 Ailsa Road Ferness Estate Ottery, Cape Town Western Cape 7800 Contact Cell No. : 082 578 3420 Contact Fax No. : 086 762 1674 Email : info@2nhealthcare.co.za

MED AIRE DUOWAVE MATTRESS SYSTEM

Item number: 820000200

The MED AIRE DuoWave alternating pressure system support offers individual cell pressure adjustment to the patient's weight thanks to its infinitely variable manual duo comfort adjustment on the unit - in both lying and sitting positions.



- Customisable thanks to the option of static soft positioning
- Stepless, manual cell pressure adjustment using a weight scale
- Bi-elastic, quilted and breathable cover
- Integrated foam pocket in the cover prevents any foam underlay used from slipping
- Quick repair of the overlay possible as air cells can be replaced individually
- Can be used for mild to high pressure ulcer risk, as well as to support pressure ulcer therapy (cat. 1-4*)

For people who can change their lying position themselves or with help on a regular basis!

Contact us today for more info.